

Self Help with Wellness Tools for parents of ASD children www.BioToolsForWellness.com

Increase Restorative Sleep

The Problem: within 3 feet of the power supply in the bedroom wall - and within 1.6 miles of a cell phone transmitter - children (and adults) have severely disturbed EEGs, decreased melatonin production and blocked activity of detoxification enzymes in the brain, gut, liver and kidneys

Diagnosis: Body Voltage meter and HF meter from BioTools

The solution: the “**Sleep Sanctuary**” technology from BioTools (evolved from the German teaching of the Institute of Baubiology).

- **Beautiful looking (and feeling) Faraday cage which decreases the EMF exposure to one tenthousands— a silver coated mosquito net and:**
- **grounding cloth**

Result: instantly normalization of the child's EEG, Heart Rate Variability and numerous hormonal, cytokine and other lab parameters. It improves blood and tissue oxygenation and tissue pH. Permanent behavioral and mental improvements in the child are dramatic and observed within days to weeks

Other helpful sleeping tools: KMT 24+ (use on feet), Valkion before sleep, toxaway foot bath at night, foot motion device

Self Help with Wellness Tools

www.BioToolsForWellness.com

Heal the brain

The problem: The astrocytes are toxic. The blood brain barrier is stuck open. The mirror neurons – responsible for normal emotional and mental responses – are toxic and dysfunctional. There is no grooming of excessive neurons, etc.

The Solutions:

The “**Valkion**”:

Swedish device which creates energized “singlet oxygen” for inhalation. This energized oxygen crosses the blood brain barrier rapidly and effectively.

The problem: in adults with neurological problems and in ASD children with Lyme disease, mycoplasma, mold and viral infections the endothelium is lined with fibrin, delaying the delivery of oxygen to the tissues from 2 seconds to 5.4 minutes. The tissues are oxygen starved. This is the most effective device to deliver oxygen and energy to the brain, to the spinal chord and the kidneys.

In our work with ASD children this is comparable with - or more effective - than HBO.

Self Help with Wellness Tools

www.BioToolsForWellness.com

The Health Light and the brain.

The visible monochromatic directional light stimulates brain function both via the visceral-cutaneous reflex and by direct penetration (up to 23 cm).

Treatment with the Health Light:

- **the deeper areas (mirror neurons, brain stem, etc.)**
- **increases bloodflow, venous and lymphatic drainage**
- **Increases metabolic activity of the brain cells**
- **activates the intrinsic detox mechanisms of the brain**
- **increases medication uptake, a most important issue that has been overlooked by most.**
- **strong reaches the cognitive areas of the brain**
- **direct anti-viral effect**
- **phototactic – it strongly attracts NK cells into the illuminated areas**

The result: the brain recovers faster, with less side effects. Appropriate treatments that did not seem to work before, can work now rapidly and profoundly.

Excellent when used together with ART. A seizure focus can be treated directly and successfully

Tools to increase the Filtration Rate of the Kidneys

- *The Problem:* the basal membranes of the kidney become inflamed and swollen as soon as mercury, lead and other toxins are filtered during chelation therapy. The nephrons stay in an allergic, poorly filtrating state, as long as the metal mobilizing medication is given
- *The solutions:*
- 1. the BioTools “**Light Shield**”. By placing the device over the kidneys normal filtrating function resumes within less then an hour. The device should be used daily during active detox and during all crises.
- 2. The Swiss “**Toxaway**” electric foot bath uses an electrochemical process to stimulate the kidney and liver meridian at the feet to effectively increase toxin elimination activity of both organs. There is also direct active lymphatic excretion of toxins across the soles of the feet. This device is superior to other comparable instruments (the secret is in the electrodes)
- *Result:* dramatic shortening of the overall treatment time and improved success rate in stubborn cases. It minimizes the feared “healing crisis” after starting a chelation program

Tools to improve phase-I and phase-II liver enzyme activity

Problem: most of our neurological clients and ASD children have impaired liver enzyme function at different locations (p450, acetylation, sulfation, glutathione-S-transferase)

- *Solution:* the BioTools “**Light Shield**”. By placing the device over the liver and kidneys
- ATP inside the cells increases significantly and the cell does its job more effectively
- turns on the cell organelles: the peroxisomes, the mitochondria, the Golgi apparatus
- enlarges the bile ducts, the arterioles and the lymph ducts in the liver by blocking the sympathetic nervous system and stimulating the parasympathetics
- Activates transport mechanisms across cell wall

Result: in spite of numerous gene polymorphisms of the detox genes, the child starts eliminating neurotoxins and heavy metals at rates comparable or superior to normal children

The KMT technology

the gentle way to handle infections in ASD children

- The instruments use 2 generators of simultaneous harmonious healthy EMF waveforms. The differential between the two frequencies creates a new rhythm and harmonics which are mathematically calculated so that they are the exact microbial inhibition frequencies published for various infectious agents. These frequencies entrain the immune system to recognize microbes which evade the system with a variety of “tricks”
- The instruments can be used like a TENS unit with electrode placement on the skin or used as a field inductor: the electrode pads are stuck to signal enhancers which are placed around the client during their sleep. It can also be used to modulate the light of the Photon Wave or the Jay Laser
- KMT 24+: use for all chronic infections as entrainment for the immune system. Fabulous results to increase effect of antimicrobial agents and medication uptake in difficult areas
- KMT 440: use for facial enhancement and anti-aging device

Lymph Massager

- This is a basic tool every health practitioner would need: increases transport of breakdown products, toxins, metabolic and cellular waste back to the subclavian vein and liver for further processing

Common uses:

- after any dental intervention (face/neck region)
- Illnesses of the eye and ear problems
- headache
- toxicity problems of brain
- regular treatment in autism for brain drainage
- Lung/liver/kidney stasis (increases drainage)
- Swollen legs/arm
- Post chemotherapy lymphatic problems
- To mobilize toxins and microbes from hidden areas

The Jae Laser

This instrument is the result of the insights of most recent biophoton physics (Prof. F.-A.Popp PhD). The green light frequency is the most effective wavelength to initiate normal action-potential and regeneration in the nervous system. The light is ideal to piggy-back other frequencies (nosodes, toxins, microbes, nutrients, homeopathics).

Most common uses:

- Brain damage
- autism
- parkinson
- MS
- insomnia
- paralysis
- infections in the spine (Bartonella, viruses)
- ALS
- Learning disorders, ADD, ADHD

Photon Wave pulsed colored light therapy

This instrument is the outcome and culmination of over 100 years of research, clinical practice and development. 2 simultaneous light sources create a holographic light beam which is modified with a 10 nanometer wide color filter. This is the most potent healing tool of energy medicine today. Treatment time varies between 6 and 10 minutes

Common uses:

- all neurological illnesses detoxification
- brain entrainment
- increase of IQ
- use to deliver medicine energetically
- chronic infections
- modification of immune system
- neurotransmitter generation and harmonization
- Insomnia
- headache
- all deep psychological interventions
- accessing the spirit world
- past life therapy

Self Help for parents of ASD children with Wellness Tools

www.BioToolsForWellness.com

- **Replicator:** this simple instrument allows to make copies of homeopathics, create auto-nosodes from urine, stool, saliva and other body fluids, create antidotes for allergens, make color remedies and more. It is based on electronic principles based on the work of the late Dr. Ludwig from Germany
- **Chi machine:** activates and normalizes sympathetic function by rhythmically pulling and tugging on the paraspinal sympathetic chain. Excellent for sleep, detoxification and body pain
- **KMT minerals:** excellent source of humic acid, minerals and magic. Best natural anti-viral (in all neurological illnesses)